# MD-Cleanser

EDTA Solution for Root Canal Cleaning and Smear Layer Removal

2020. 10. 16. DM-51-16

Manufacturer

**META BIOMED CO., LTD.**

270, Osongsangmyeong 1-ro,  
Osong-eup, Heungdeok-gu, Cheongju-si,  
Chungcheongbuk-do, Republic of Korea

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Fax +49-208-30991999

en INSTRUCTION FOR USE

**EDTA Solution for Root Canal Cleaning and Smear Layer Removal**

**Excellent Smear Layer Removal**

Presence of calcifications in the canal, too narrow canals in posterior teeth, removal of organic in root canal, and smear layer removal.

Application of MD - Cleanser For Biomechanical Preparation of the canal Apply MD - Cleanser by depositing a few drops into the root canal using a syringe, and then carefully pump the solution into the canal with a fine root canal instrument.

Continue instrumentation of the canal while being soaked in the solution.

Re - apply MD - Cleanser when necessary until canal shaping and cleaning is complete.

Remove MD - Cleanser from the canal by irrigation. For Smear Layer Removal Irrigate the canal with MD - Cleanser, followed immediately by thorough irrigation with Water to rinse off MD - Cleanser. Aspirate excess fluid and dry the canal.

**Properties**

MD - Cleanser functions by forming a calcium - chelate solution with the calcium ions of dentin. It makes the dentin brittle and easier to instrument, and produces the cleanest dentinal walls after smear layer

**Special precautions for use**

Use MD - Cleanser with discretion to prevent perforation. Avoid extrusion of the solution into the apex. Keep tightly capped at all times. For dental use only. Keep out of children's reach.

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